



FAST RESULT It takes only three minutes to change the color of the test from yellow to blue.

NO INCUBATION Unlike most of the rapid urease tests, AMA® RUT 10 detects the urease that is present in the biopsy specimen at the time it is taken, not the urease that can be produced hours later. Thus, saving time on incubation, AMA® RUT 10 also provides better specificity compared to most of the gel-based urease tests.

ACCURACY APPROVED AMA® RUT 10 is an accurate point-of-care test: it proved to be highly sensitive (93%) and specific (95%).

CE MARKED AMA® RUT meets the requirements of the In-vitro Diagnostics Directive 98/79/EC and is officially registered in European Community.

USE THE SPECIMEN FURTHER The same biopsy specimen can be examined with other methods such as culture, histology etc. straight after the test.

EASY TO STORE Shelf-life of the test is 24 months. No special storage conditions are required.

EVIDENT COLOR CHANGE While other urease tests change color from yellow to red, AMA® RUT 10 reacts positive by changing its color to blue. This is easier to notice, because the color change differs from the color of the biopsy specimen.



AMA RAPID UREASE TEST AMA® RUT 10

The Up-to-date In Vitro Diagnostics




AMA Co Ltd
199034, 4-6, 17 liniya V.O.,
Saint-Petersburg, Russia
tel./fax 007-812-380-7699,
tel. 007-812-321-7501
www.amamed.ru; ama@sp.ru

AMA® Rapid Urease Test 10 (AMA® RUT 10) is the new format of the ready-to-use test for detection of *H. pylori* urease activity in biopsy specimen obtained by gastroscopy from patients with suspected *H. pylori* infection.

ECONOMY




Multiple tests on a single AMA® RUT 10 device. Help to reduce the price per test.




AMA® RUT 10 with unconsumed parts can be used further

As you open one segment, other segments remain hermetically sealed and can be further used or stored.




*200 pcs=2000 tests

Compact slim design of the new AMA® RUT 10



Saves doctor's time

The same biopsy specimen can be used for further diagnostics by other methods such as culture, histology etc. after AMA® RUT 10 has been applied. Thus, thanks to AMA® RUT 10 doctor saves time on collecting the samples.



FOR FREE SAMPLES

